

DX-G

NEXT-GENERATION CR SYSTEM

The DX-G digitizer unites superb image quality with a drop-and-go buffer-based workflow and enables a potential reduction in patient dose. It offers the unprecedented convenience of being able to combine standard phosphor plates and needle-based detectors.

- State-of-the-art image quality, with potential dose reduction
- Drop-and-go cassette buffer
- Broad range of applications
- Both needle-based detectors and standard phosphor plates

The next-generation in CR for general radiography departments, the DX-G digitizer unites superb image quality with the convenience of supporting both standard phosphor plates and needle-based detectors. The exclusive DirectriX detector technology offers the potential for a significant patient dose reduction. With a user-friendly drop-and-go buffer that can handle a mix of five cassettes of different sizes, workflow is smoother and more productive. The DX-G can be used as a centralized or decentralized digitizer in the radiography department, supporting a broad range of applications. In a centralized environment, it can serve multiple rooms. At the same time, its small footprint means it can be placed in any available space.


State-of-the art image quality, with potential dose reduction

By supporting both standard phosphor plates and needle-based detectors, the DX-G unites complete convenience with high image quality, while leveraging the radiography department's existing investments. With standard phosphor plates, the DX-G delivers excellent image quality. When used with DirectriX needle-based detectors, however, the DX-G provides superb image quality with a much higher Detective Quantum Efficiency (DQE). This state-of-the-art image quality offers the potential to reduce patient dose.





Needle-based
detector



Powder
phosphor plate

Broad range of applications

The combination of needle-based detectors, standard phosphor plates with specific cassettes and image resolution mode make the DX-G ideal for a broad range of applications:

- General radiography
- Orthopedics - extremities
- Dental
- Pediatrics and neonatal
- Full Leg / Full Spine

It offers two different image resolution modes:
100 µm pixel pitch (10 pixels/mm) and 150 µm pixel pitch (6.7 pixels/mm).

Optimal productivity and smooth workflow

The drop-and-go buffer and fast preview eliminate waiting times and facilitate a continuous workflow within the department.

The five-cassette drop-and-go buffer can handle a mix of different sizes of both needle-based detectors and standard phosphor plates. The automatic cassette handling makes DX-G highly productive and user-friendly.

Using DX-G as a central digitizer in the radiography department, multiple examination rooms can be supported. With its small footprint, it can fit into the tightest spaces, including the X-ray room or even a narrow corridor.

The right choice

To eliminate any confusion, needle-based detector cassettes are gray, while standard phosphor plate cassettes are orange, so that there is little chance of the user making a mistake when selecting the desired cassette. Each plate has an embedded memory that stores the data entered during identification by no-touch radiofrequency tagging. Thus, the identification data and images are linked from the beginning throughout the entire digital processing system.

SAFETY

Region	Safety	EMC	Laser
EUROPE	IEC 60601-1:1988 + A1:1991: + A2:1995	EN 60601-1-2:2007 EN 300 330 2 V1.1.1:2001 EN 301 489 V1.3.1:2001	60825-1:1993 + A1:1997 + A2:2001
USA	UL60601-1:2003	FCC part 15	CFR parts 1040.10 and 1040.11
CANADA	CSA C 22.2 No.601.1: 1990 + S1:1994 + A2:1998	CSA C 22.2 No. 601.1.2	CSA-E60825-1-03

Cassettes for
needle-based
detectorsCassettes for
standard phosphor
plates

DETECTORS

Needle-based detector

- CR HD5.0 General SR
- CR HD5.0 General

Size

35 x 43 cm (13.77 x 16.92")
 35 x 43 cm (13.77 x 16.92")
 24 x 30 cm (9.44 x 11.81")
 18 x 24 cm (7.08 x 9.44")
 15 x 30 cm (5.90 x 11.81")
 35 x 43 cm (13.77 x 16.92")
 24 x 30 cm (9.44 x 11.81")
 18 x 24 cm (7.08 x 9.44")
 35 x 43 cm (13.77 x 16.92")

Spatial resolution

6.7 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm

Pixel matrix

2272 x 2800
 3408 x 4200
 2256 x 2880
 1656 x 2280
 1344 x 2880
 3408 x 4200
 2256 x 2880
 1656 x 2280
 3408 x 4368

Standard phosphor plate

- CR MD4.0R General SR
- CR MD4.0R General
- CR MD4.0R FLFS SR

Size

35 x 43 cm (13.77 x 16.92")
 35 x 35 cm (13.77 x 13.77")
 35 x 43 cm (13.77 x 16.92")
 35 x 35 cm (13.77 x 13.77")
 24 x 30 cm (9.44 x 11.81")
 18 x 24 cm (7.08 x 9.44")
 15 x 30 cm (5.90 x 11.81")
 35 x 43 cm (13.77 x 16.92")

Spatial resolution

6.7 pixels/mm
 6.7 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm
 10 pixels/mm

Pixel matrix

2320 x 2832
 2320 x 2320
 3480 x 4248
 3480 x 3480
 2328 x 2928
 1728 x 2328
 1440 x 2928
 3480 x 4392

Technical Specifications

GENERAL

Drop-and-go cassette buffer

- 5 cassettes of mixed sizes input buffer and 5 cassettes of mixed sizes output buffer

Throughput

- 35 x 43 cm (14 x 17") = approx. 83 plates/h

Display for status and error indication

- LCD touchscreen
- LED status indicator

Greyscale resolution

- Output to processor:
16 bits/pixel square root compressed

Dimensions and weight

- Covered floor space (W x D x H):
660 x 510 x 1230 mm (26 x 20 x 48.4")
- Output buffer included (W x D x H):
1150 x 510 x 1230 mm (43.5 x 20 x 48.4")
- Weight: approx. 180 kg (397 lb)

Configuration requirements

- NX
- ID tablet
- CR HD5.0 Detectors and Cassettes
- CR MD4.0R Plates and Cassettes

Power

- 220 - 240 V/50-60 Hz
Standby 87 W, peak 590 W, fuse 16 A
- 120 V/60 Hz (USA)
Standby 92 W, peak 621 W, fuse 15 A
- 100 V/60 Hz (Japan)
Standby 92 W, peak 621 W, fuse 15 A

Environmental Requirements

DX-G digitizer

- Temperature: +15 ~ +30° C (+59 ~ +86° F)
- Humidity: 15 ~ 75% Rh
- EMC compliant with IEC 60601-1-2
- Rate of change of temperature:
0.5° C/minute (0.9° F)

Transport details

- Temperature: -25 ~ +55° C (-13 to 131° F),
-25° C for max. 72 hours, +55° C for max. 96 hours
- Humidity: 5 ~ 95% Rh

Environmental effects

- Noise level: max. 65 dB (A)
- Heat dissipation: standby 92 W,
continuous operation 242 W

SAFETY

Approvals

- ETL classified CUS, CE

For more information on Agfa, please visit our website on www.agfa.com ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. DirectriX and the DirectriX logo are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2018 Agfa NV
All rights reserved
Published by Agfa NV
Septestraat 27 - 2640 Mortsel
Belgium

50C6K US 00201806