

AGFA HEALTHCARE DICOM Conformance Statement



IMPAX HeartStation 1.7

Document No. 001494
Revision 1.0
Livelihood NodeID: 48739041

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp
--	---	---

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel
Belgium

tel: 32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment by Agfa personnel.

Copyright © March, 15
Agfa HealthCare
All rights reserved

Conformance Statement Overview

IMPAX HeartStation will create and consume DICOM Waveform and DICOM Encapsulated PDF objects. The following services are supported by HeartStation:

Table 1.1-1: Network Services Supported

DICOM SOP Class Name	User of Service (SCU)	Provider of Service (SCP)
Transfer		
12-lead ECG Waveform Storage	Yes	Yes
Encapsulated PDF Storage	Yes	Yes
Workflow Management		
Modality Worklist Information Model - FIND	Yes	No

Table of Contents

Document Information	2
Conformance Statement Overview	3
Table of Contents	4
1 Introduction	5
1.1 Revision Record	5
1.2 Purpose and Intended Audience of this Document	5
1.3 General Remarks	5
1.3.1 Integration and Validation Activities	5
1.3.2 Future Evolution	5
1.4 Acronyms and Abbreviations	6
1.5 Related Documents	6
2 Networking	7
2.1 Implementation Model	7
2.1.1 Application Data Flow Diagram	7
2.1.2 Functional Definitions of AE's	8
2.1.2.1 Functional Capability of SCP	8
2.1.2.2 Functional Capability of SCU	8
2.2 AE Specifications	8
2.2.1 HeartStation Specification	8
2.2.1.1 SOP Classes Supported	8
2.2.1.2 Association Establishment Policies	8
2.2.1.2.1 General	8
2.2.1.2.2 Number of Associations	8
2.2.1.2.3 Asynchronous Nature	9
2.2.1.2.4 Implementation Identifying Information	9
2.2.1.3 Association Initiation Policies	9
2.2.1.3.1 Activity – SCU	9
2.2.1.3.1.1 Description and Sequencing of Activity	9
2.2.1.3.1.2 Proposed Presentation Contexts	9
2.2.1.3.2 SOP Specific Conformance for Encapsulated PDF IOD Storage	10
2.2.1.3.3 SOP Specific Conformance for Modality Worklist Queries	11
2.2.1.3.4 SOP Specific Conformance for 12-lead ECG Waveform Storage IOD Storage	12
2.2.1.3.5 Activity – SCP	18
2.2.1.3.5.1 Description and Sequencing of Activity	18
2.2.1.3.6 Proposed Presentation Contexts	18
2.2.1.3.7 SOP Specific Conformance for Encapsulated PDF IOD Storage	19
2.2.1.3.8 SOP Specific Conformance for 12-lead ECG Waveform Storage	19
2.3 Network Interfaces	21
2.3.1 Physical Medium Support	21
2.4 Configuration	22
2.4.1 Configuration Parameters	22
3 Support for Extended Character Sets	23
4 Security	24
4.1 Security Profiles	24
4.2 Association Level Security	24

4.3	Application Level Security	24
-----	----------------------------------	----

1 INTRODUCTION

1.1 Revision Record

DICOM Conformance Statement IMPAX HeartStation 1.7		
Revision Number	Date	Reason for Change
1.0	February 13, 2015	Initial version for HeartStation 1.7

1.2 Purpose and Intended Audience of this Document

This document is a DICOM Conformance Statement for the DICOM Services of the HeartStation product. It is written according to part PS 3.2 of Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-3.16, 2008.

The user of this document is involved with system integration and/or software design. We assume that the reader is familiar with the terminology and concepts that are used in the DICOM 3.0 standard and the IHE Technical Framework.

Readers not familiar with DICOM 3.0 terminology should first read the appropriate parts of the DICOM standard itself, prior to reading this conformance statement.

1.3 General Remarks

1.3.1 Integration and Validation Activities

The integration of any device into a system of interconnected devices goes beyond the scope of the DICOM 3.0 standard and this conformance statement when *interoperability* is desired. The responsibility for analyzing the applications requirements and developing a solution that integrates the Agfa equipment with other vendors' systems is the user's responsibility and should not be underestimated.

In some circumstances it might be necessary to perform a validation to make sure that functional interoperability between the Agfa equipment and non-Agfa devices works as expected. The user should ensure that any non-Agfa provider accepts responsibility for any validation required for their connection with the Agfa equipment.

1.3.2 Future Evolution

As the DICOM 3.0 standard evolves to meet the user's growing requirements and to incorporate new features and technologies, Agfa will follow the evolution of the standard. This evolution of the standard may require changes to devices that have implemented DICOM 3.0. The user should ensure that any non-Agfa provider, who connects with Agfa devices, also plans for future evolution of the DICOM standard. A refusal to do so may result in the loss of functionality and/or connectivity between the different products.

1.4 Acronyms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard. Abbreviations and terms are as follows:

AE	DICOM Application Entity
DICOM	Digital Imaging and Communications in Medicine
GUI	Graphical User Interface
IHE	Integrating the Healthcare Enterprise
IOD	(DICOM) Information Object Definition
PACS	Picture Archive and Communications System
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM server)
SOP	DICOM Service-Object Pair
UID	Unique Identifier
UTF-8	Unicode Transformation Format - 8
VR	Value Representation

1.5 Related Documents

- ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) V3.0.2008
- IHE Radiology Technical Framework Revision 5.5 – Final Text, November 20, 2003
- Health Level Seven standard version 2.3.1

2 NETWORKING

2.1 Implementation Model

2.1.1 Application Data Flow Diagram

The Application Data Flow Diagram in Figure 2.1-1 depicts the DICOM data flow to and from HeartStation. The tail of the arrow between HeartStation and the remote real world activity indicates the party that initiates the association negotiation.




Figure 2.1-1: Functional Overview – Application Data Flow

2.1.2 Functional Definitions of AE's

This section will provide a brief description of the various application entities that are included as a part of HeartStation.

2.1.2.1 Functional Capability of SCP

The HeartStation SCP is capable of receiving 12-lead ECG Waveform Storage and Encapsulated PDF Storage requests.

2.1.2.2 Functional Capability of SCU

The HeartStation SCU is capable of sending Encapsulated PDF Storage and 12-lead ECG Waveform Storage messages. The HeartStation SCU is capable of sending Modality Worklist Queries.

2.2 AE Specifications

2.2.1 HeartStation Specification

2.2.1.1 SOP Classes Supported

This Application Entity provides Standard Conformance to the following SOP Classes:

Table 2.2-1: SOP Classes for HeartStation

SOP Class Name	SOP Class UID	SCU	SCP
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Yes	Yes
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Yes	Yes
Modality Worklist Information Model - FIND	1.2.840.10008.5.1.4.31	Yes	No

2.2.1.2 Association Establishment Policies

2.2.1.2.1 General

The DICOM standard Application context shall be specified.

Table 2.2-2: DICOM Application Context

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

2.2.1.2.2 Number of Associations

The maximum number of simultaneous associations accepted by HeartStation is not constrained. There is no inherent limit to the number of associations other than limits imposed by the computer operating system.

2.2.1.2.3 Asynchronous Nature

HeartStation allows a single outstanding operation on any association. Therefore, HeartStation does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

2.2.1.2.4 Implementation Identifying Information

HeartStation will respond with the following implementation identifying parameters:

Table 2.2-3: DICOM implementation Class and Version for Verification

SCP Identification Parameters	
Implementation Class UID	1.2.40.0.13.1.1
Implementation Version Name	dcm4che-1.4.28
SCU Identification Parameters	
Implementation Class UID	1.2.40.0.13.1.1
Implementation Version Name	dcm4che-1.4.28

2.2.1.3 Association Initiation Policies

2.2.1.3.1 Activity – SCU

2.2.1.3.1.1 Description and Sequencing of Activity

HeartStation will act as an SCU to store an Encapsulated PDF Storage objects and 12-lead ECG Waveform Storage objects. HeartStation will act as an SCU to perform Modality Worklist queries.

2.2.1.3.1.2 Proposed Presentation Contexts

Table 2.2-4: Presentation Contexts Proposed

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None
Modality Worklist Information Model - FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

2.2.1.3.2 SOP Specific Conformance for Encapsulated PDF IOD Storage

HeartStation creates an Encapsulated PDF object with the following DICOM attributes. All attributes conform to the Encapsulated PDF IOD class specification.

Table 2.2-5: Encapsulated PDF IOD Storage Attributes

Module	Attribute Name	Tag	Value
Patient	Patient's Name	(0010,0010)	ECG patient first and last names
	Patient ID	(0010,0020)	ECG patient id
	Patient's Birth Date	(0010,0030)	ECG birth date
	Patient's Sex	(0010,0040)	ECG patient gender
General Study	Study Instance UID	(0020,000D)	Generated by HeartStation using the following prefix : 2.25
	Study Date	(0008,0020)	ECG acquisition date
	Study Time	(0008,0030)	ECG acquisition time
	Referring Physician's Name	(0008,0090)	From ECG or unvalued
	Study ID	(0020,0010)	From ECG, or unvalued
	Accession Number	(0008,0050)	From HIS, or unvalued
Encapsulated Document Series	Modality	(0008,0060)	"ECG"
	Series Instance UID	(0020,000E)	Generated by HeartStation using the following prefix : 2.25
	Series Number	(0020,0011)	"1"
	Series Description	(0008,103E)	"ECG Report - Embedded PDF"
General Equipment	Manufacturer	(0008,0070)	"Agfa"
SC Equipment	Conversion Type	(0008,0064)	"SD"
Encapsulated Document	Instance Number	(0020,0013)	"1"
	Content Date	(0008,0023)	Same value as Study Date
	Content Time	(0008,0033)	Same value as Study Time
	Acquisition DateTime	(0008,002A)	From ECG
	Burned In Annotation	(0028,0301)	"YES"
	Document Title	(0042,0010)	"ECG Report"
	> Concept Code Name Sequence	(0040,A043)	-
	>> Code Value	(0008,0100)	"ECG Report"
	>> Coding Scheme Designator	(0008,0102)	"ECG Report"
	>> Code Meaning	(0008,0104)	"1.1"
	>> Code Scheme Version	(0008,0103)	"1"
	MIME Type of Encapsulated Document	(0042,0012)	"application/pdf"
	Encapsulated	(0042,0011)	Encapsulated Document
	Verification Flag	(0040,A493)	If ECG was confirmed, then value is VERIFIED, if ECG was not confirmed value is "UNVERIFIED"
SOP Common	SOP Class UID	(0008,0016)	"1.2.840.10008.5.1.4.1.1.104.1"
	SOP Instance UID	(0008,0018)	Generated by HeartStation using the following components: AGFA prefix: 1.2.124.113532
	Instance Creation Date	(0008,0012)	Date the SOP instance was created

Module	Attribute Name	Tag	Value
	Instance Creation Time	(0008,0013)	Time the SOP instance was created
	Specific Character Set	(0008,0005)	"ISO_IR 100"
	Predecessor Document Sequence	(0040,A360)	
	> Study Instance UID	(0020,000D)	Study Instance UID of the previous object
	> Referenced Series Sequence	(0008,1115)	
	>> Series Instance UID	(0020,000E)	Series Instance UID of the previous object
	>> Referenced SOP Sequence	(0008,1199)	Sequence containing an item for each previous SOP instance
	>>> Referenced SOP Instance UID	(0008,1155)	SOP Instance UID of the previous object
	>>> Referenced SOP Class UID	(0008,1150)	SOP Class UID of the previous object
	>>> Purpose of Reference Sequence	(0040,A170)	
	>>>> Code Value	(0008,0100)	"121361"
	>>>> Coding Scheme Designator	(0008,0102)	"DCM"
	>>>> Coding Scheme Designator	(0008,0104)	"Addended Report"

2.2.1.3.3 SOP Specific Conformance for Modality Worklist Queries

HeartStation provides standard conformance to the DICOM Modality Worklist Information Model Class as an SCU. The DICOM attributes that are utilized by HeartStation when issuing a DICOM Modality Worklist request are listed in Table 2.2-6. Queries are only by Patient ID (0010,0020)

Table 2.2-6: Attributes for the Modality Worklist Information Model

Attribute	Tag
Accession Number	(0008,0050)
Admission ID	(0038,0010)
Current Patient Location	(0038,0300)
Patient's Birth Date	(0010,0030)
Patient's Name	(0010,0010)
Patient's Sex	(0010,0040)
Patient's Race	(0010,2160)
Patient's State	(0038,0500)
Patient ID	(0010,0020)
Requested Procedure Description	(0032,1060)
Requested Procedure Code	(0032,1064)
Requesting Physician	(0032,1032)
Scheduled Procedure Step Start Date	(0040,0002)
Scheduled Procedure Step Start Time	(0040,0003)
Scheduled Procedure Step Status	(0040,0020)
Placer Order Number	(0040,2016)
Filler Order Number	(0040,2017)

2.2.1.3.4 SOP Specific Conformance for 12-lead ECG Waveform Storage IOD Storage

HeartStation creates an 12-lead ECG Waveform object with the following DICOM attributes. All attributes conform to the 12-lead ECG Waveform IOD class specification.

Table 2.2-7: 12-lead ECG Waveform Storage Storage Attributes

Module	Attribute Name	Tag	Value
Patient	Patient's Name	(0010,0010)	ECG patient first and last names
	Patient ID	(0010,0020)	ECG patient id
	Patient's Birth Date	(0010,0030)	ECG birth date
	Patient's Sex	(0010,0040)	ECG patient gender
General Study	Study Instance UID	(0020,000D)	Generated by HeartStation with the following prefix: 2.25
	Study Date	(0008,0020)	ECG acquisition date
	Study Time	(0008,0030)	ECG acquisition time
	Referring Physician's Name	(0008,0090)	From ECG or unvalued
	Study ID	(0020,0010)	From ECG, or unvalued
	Accession Number	(0008,0050)	From HIS, or unvalued
	Study Description	(0008,1030)	From ECG, or unvalued
	Patient Size	(0010,1020)	From ECG, or unvalued
	Patient Weight	(0010,1030)	From ECG, or unvalued
General Series	Modality	(0008,0060)	"ECG"
	Series Instance UID	(0020,000E)	Generated by HeartStation with the following prefix: 2.25
	Series Number	(0020,0011)	"1"
	Series Description	(0008,103E)	"ECG Waveform "
	Operators Name	(0008,1070)	From ECG, or unvalued
General Equipment	Manufacturer	(0008,0070)	From ECG, or unvalued
	Manufacturer Model Name	(0008,1090)	From ECG, or unvalued
	Station Name	(0008,1010)	From ECG, or unvalued
	Institution Name	(0008,0080)	From ECG, or unvalued
	Institutional Department Name	(0008,1040)	From ECG, or unvalued
Waveform Identification	Instance Number	(0020,0013)	"1"
	Content Date	(0008,0023)	Same value as Study Date
	Content Time	(0008,0033)	Same value as Study Time
	Acquisition DateTime	(0008,002A)	Same value as Study Date/Time
Waveform	Waveform Sequence	(5400,0100)	First item will contain the original waveform group, and the second item will contain the median beat group.
	> Multiplex Group Time Offset	(0018,1068)	"0"
	> Trigger Time Offset	(0018,1069)	"0"
	> Waveform Originality	(003A,0004)	"ORIGINAL" for rhythm group, "DERIVED" for median beat group
	> Number of Waveform Channels	(003A,0005)	Number is leads in this group. Either "12" or "15"

Module	Attribute Name	Tag	Value
	> Number of Waveform Samples	(003A,0010)	Number of samples in each lead
	> Sampling Frequency	(003A,001A)	Original sample rate from ECG. Either "250" or "500"
	> Multiplex Group Label	(003A,0020)	"RHYTHM" for rhythm group, or "MEDIAN BEAT" for median beat group
	> Channel Definition Sequence	(003A,0200)	Sequence of items, with one item per channel
	>> Channel Source Sequence	(003A,0208)	
	>>> Code Value	(0008,0100)	The code value for the lead, one of: Lead I : "5.6.3-9-1", Lead II: "5.6.3-9-2", Lead III: "5.6.3-9-61", Lead AVR: "5.6.3-9-62", Lead AVL: "5.6.3-9-63", Lead AVF: "5.6.3-9-64", Lead V1: "5.6.3-9-3", Lead V2: "5.6.3-9-4", Lead V3: "5.6.3-9-5", Lead V4: "5.6.3-9-6", Lead V5: "5.6.3-9-7", Lead V6: "5.6.3-9-8", Lead V7: "5.6.3-9-9", Lead V8: "5.6.3-9-66", Lead V9: "5.6.3-9-67", Lead V3R: "5.6.3-9-11", Lead V4R: "5.6.3-9-12", Lead E1: "5.6.3-9-75", Lead E2: "5.6.3-9-76", Lead E3: "5.6.3-9-77"
	>>> Coding Scheme Designator	(0008,0102)	"SCPECG"
	>>> Coding Scheme Version	(0008,0103)	"1.3"

Module	Attribute Name	Tag	Value
	>>> Code Meaning	(0008,0104)	The code meaning for the lead, one of: Lead I: "Lead I (Einthoven)", Lead II: "Lead II", Lead III: "Lead III", Lead aVR: "Lead aVR", Lead aVL: "Lead aVL", Lead aVF: "Lead aVF", Lead V1: "Lead V1", Lead V2: "Lead V2", Lead V3: "Lead V3", Lead V4: "Lead V4", Lead V5: "Lead V5", Lead V6: "Lead V6", Lead V7: "Lead V7", Lead V8: "Lead V8", Lead V9: "Lead V9", Lead V3R: "Lead V3R", Lead V4R: "Lead V4R", Lead E1: "Lead E1", Lead E2: "Lead E2", Lead E3: "Lead E3"
	>> Channel Sensitivity	(003A,0210)	Nominal numeric value of a unit quality of a sample. This is the LSB in μV . Will be one of: 1, 2.44, 2.5, 4.88, 5, 1.22, 1.25, 2, 4, 8, 9.76, 10, 20
	>> Channel Sensitivity Units	(003A,0211)	A coded descriptor of the Units of measure for the channel sensitivity
	>>> Code Value	(0008,0100)	" μV "
	>>> Coding Scheme Designator	(0008,0102)	"UCUM"
	>>> Coding Scheme Version	(0008,0103)	"1.4"
	>>> Code Meaning	(0008,0104)	microvolt
	>> Channel Sensitivity Correction Factor	(003A,0212)	"1"
	>> Channel Baseline	(003A,0213)	"0"
	>> Channel Sample Skew	(003A,0215)	"0"
	>> Waveform Bits Stored	(003A,021A)	"16"
	>> Filter Low Frequency	(003A,0220)	One of: 0.05, 0.15, or 0.50
	>> Filter High Frequency	(003A,0221)	One of: 40, 100, or 150
	>> Notch Filter Frequency	(003A,0222)	Absent if no notch filter applied otherwise one of: 50 or 60
	> Waveform Bits Allocated	(5400,1004)	"16"
	> Waveform Sample Interpretation	(5400,1006)	"SS"
	> Waveform Data	(5400,1010)	Encoded ECG sample data
Acquisition Context	Acquisition Context Sequence	(0040,0555)	
	> Value Type	(0040,A040)	"CODE"
	> Concept Name Code Sequence	(0040,A043)	
	>> Code Value	(0008,0100)	"5.4.5-33-1"

Module	Attribute Name	Tag	Value
	>> Coding Scheme Designator	(0008,0102)	"SCPECG"
	>> Coding Scheme Version	(0008,0103)	"1.3"
	>> Code Meaning	(0008,0104)	"Electrode Placement"
	> Concept Code Sequence	(0040,A168)	
	>> Code Value	(0008,0100)	"5.4.5-33-1-1"
	>> Coding Scheme Designator	(0008,0102)	"SCPECG"
	>> Coding Scheme Version	(0008,0103)	"1.3"
	>> Code Meaning	(0008,0104)	"Standard 12-lead positions: limb leads placed at extremities"
	The following is present if Systolic Blood Pressure is available in ECG		
	> Concept Name Code Sequence	(0040,A043)	
	>> Code Value	(0008,0100)	"F-008EC"
	>> Coding Scheme Designator	(0008,0102)	"SRT"
	>> Code Meaning	(0008,0104)	"Systolic Blood Pressure"
	> Measurement Units Code Sequence	(0040,08EA)	
	>> Code Value	(0008,0100)	"mmHg"
	>> Coding Scheme Designator	(0008,0102)	"UCUM"
	>> Coding Scheme Version	(0008,0103)	"1.4"
	>> Code Meaning	(0008,0104)	"millimeters Hg"
	> Value Type	(0040,A040)	"NUM"
	> Numeric Value	(0040,A30A)	Systolic Blood Pressure Value
	The following is present if Diastolic Blood Pressure is available in ECG		
	> Concept Name Code Sequence	(0040,A043)	
	>> Code Value	(0008,0100)	"F-008ED"
	>> Coding Scheme Designator	(0008,0102)	"SRT"
	>> Code Meaning	(0008,0104)	"Diastolic Blood Pressure"
	> Measurement Units Code Sequence	(0040,08EA)	
	>> Code Value	(0008,0100)	"mmHg"
	>> Coding Scheme Designator	(0008,0102)	"UCUM"
	>> Coding Scheme Version	(0008,0103)	"1.4"
	>> Code Meaning	(0008,0104)	"millimeters Hg"
	> Value Type	(0040,A040)	"NUM"
	> Numeric Value	(0040,A30A)	Diastolic Blood Pressure Value
Waveform Annotation	Waveform Annotation Sequence	(0040,B020)	Sequence of annotation items
	The following repeats for each ECG interpretation statement		
	> Annotation Group Number	(0040,A180)	"0"
	> Unformatted Text Value	(0070,0006)	The ECG interpretation statement
	> Referenced Waveform Channels	(0040,A0B0)	"1\0"

Module	Attribute Name	Tag	Value
	The following repeats for each measurement		
	> Concept Name Code Sequence	(0040,A043)	
	>> Code Value	(0008,0100)	One of the following depending on the type of measurement: PR Interval: "5.13.5-7", QRS Duration: "5.13.5-9", QT Interval: "5.13.5-11", QTc Interval: "5.10.2.5-5", P Axis: "5.10.3-11", QRS Axis: "5.10.3-13", T Axis: "5.10.3-15", Ventricular Heart Rate: "5.10.2.5-1"
	>> Coding Scheme Designator	(0008,0102)	"SCPECG"
	>> Coding Scheme Version	(0008,0103)	"1.3"
	>> Code Meaning	(0008,0104)	One of the following depending on the type of measurement: PR Interval: "PR Interval", QRS Duration: "QRS Duration", QT Interval: "QT Interval", QTc Interval: "QTc Interval", P Axis: "P Axis", QRS Axis: "QRS Axis", T Axis: "T Axis", Ventricular Heart Rate: "Ventricular Heart Rate"
	> Measurement Units Code Sequence	(0040,08EA)	
	>> Code Value	(0008,0100)	Depending on the measurement, will be one of: PR Interval: "ms", QRS Duration: "ms", QT Interval: "ms", QTc Interval: "ms", P Axis: "deg", QRS Axis: "deg", T Axis: "deg", Ventricular Heart Rate: "{H.B.}/min"
	>> Coding Scheme Designator	(0008,0102)	UCUM
	>> Coding Scheme Version	(0008,0103)	1.4
	>> Code Meaning	(0008,0104)	Depending on the measurement, will be one of: PR Interval: "millisecond", QRS Duration: "millisecond", QT Interval: "millisecond", QTc Interval: "millisecond", P Axis: "degrees", QRS Axis: "degrees", T Axis: "degrees", Ventricular Heart Rate: "heart beats per minute"

Module	Attribute Name	Tag	Value
	> Referenced Waveform Channels	(0040,A0B0)	"1\0"
	> Annotation Group Number	(0040,A180)	"1"
	> Numeric Value	(0040,A30A)	The measured value
	The follow repeats for each pacing marker (if present)		
	> Concept Name Code Sequence	(0040,A043)	
	>> Code Value	(0008,0100)	5.10.1.2
	>> Coding Scheme Designator	(0008,0102)	"SCPECG"
	>> Coding Scheme Version	(0008,0103)	"1.3"
	>> Code Meaning	(0008,0104)	Pacemaker spike
	> Referenced Sample Positions	(0040,A132)	Sample position of pacemaker spike
	> Referenced Waveform Channels	(0040,A0B0)	"1\0"
	> Annotation Group Number	(0040,A180)	"3"
	> Numeric Value	(0040,A30A)	The measured value
SOP Common	SOP Class UID	(0008,0016)	"1.2.840.10008.5.1.4.1.1.9.1.1"
	SOP Instance UID	(0008,0018)	Generated by HeartStation with the following prefix: 2.25
	Instance Creation Date	(0008,0012)	Date the SOP instance was created
	Instance Creation Time	(0008,0013)	Time the SOP instance was created
	Specific Character Set	(0008,0005)	"ISO_IR 100"
	Predecessor Document Sequence	(0040,A360)	
	> Study Instance UID	(0020,000D)	Study Instance UID of the previous object
	> Referenced Series Sequence	(0008,1115)	
	>> Series Instance UID	(0020,000E)	Series Instance UID of the previous object
	>> Referenced SOP Sequence	(0008,1199)	Sequence containing an item for each previous SOP instance
	>>> Referenced SOP Instance UID	(0008,1155)	SOP Instance UID of the previous object
	>>> Referenced SOP Class UID	(0008,1150)	SOP Class UID of the previous object
	>>> Purpose of Reference Sequence	(0040,A170)	
	>>>> Code Value	(0008,0100)	"121361"
	>>>> Coding Scheme Designator	(0008,0102)	"DCM"
	>>>> Coding Scheme Designator	(0008,0104)	"Addended Report"
Private Tags		(2633,0010)	"Agfa IMPAX HeartStation Waveform"
		(2633,1000)	HeartStation user text field #1 label
		(2633,1001)	HeartStation user text field #1 value
		(2633,1002)	HeartStation user text field #2 label
		(2633,1003)	HeartStation user text field #2 value

Module	Attribute Name	Tag	Value
		(2633,1004)	HeartStation user text field #3 label
		(2633,1005)	HeartStation user text field #3 value
		(2633,1006)	HeartStation user text field #4 label
		(2633,1007)	HeartStation user text field #4 value
		(2633,1010)	HeartStation account number
		(2633,1011)	HeartStation order number
		(2633,1012)	HeartStation ECG severity code name, if present
		(2633,1020)	HeartStation ECG department name
		(2633,1021)	HeartStation ECG room number
		(2633,1030)	HeartStation ECG DX value #1, if present
	
		(2633,1034)	HeartStation ECG DX value #5, if present
		(2633,1030)	HeartStation ECG RX value #1, if present
	
		(2633,1034)	HeartStation ECG RX value #5, if present
		(2633,1050)	HeartStation current ECG study state
		(2633,1051)	HeartStation current ECG study state activation date time
		(2633,1052)	HeartStation current ECG study state applied by name, if present
		(2633,1053)	HeartStation current ECG study state signing proxy name, if present
		(2633,1055)	HeartStation current ECG status statement
		(2633,1060)	HeartStation ECG trace sensitivity display in mm per mV. One of: "5," "10", "20", "40"

2.2.1.3.5 Activity – SCP

2.2.1.3.5.1 Description and Sequencing of Activity

HeartStation will act as an SCP to receive 12-lead ECG Waveform Storage and Encapsulated PDF Storage objects.

2.2.1.3.6 Proposed Presentation Contexts

Table 2.2-8: Presentation Contexts Proposed

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCP	None
--------------------------	-------------------------------	--	--	-----	------

2.2.1.3.7 SOP Specific Conformance for Encapsulated PDF IOD Storage

HeartStation accepts an Encapsulated PDF object and uses the following DICOM attributes. All attributes conform to the Encapsulated PDF IOD class specification.

Table 2.2-9: Encapsulated PDF IOD Storage Attributes

Module	Attribute Name	Tag	Notes
Patient	Patient's Name	(0010,0010)	
	Patient ID	(0010,0020)	
	Patient's Birth Date	(0010,0030)	Optional
	Patient's Sex	(0010,0040)	Optional
Patient Demographics	Patient's Height	(0010,1020)	Optional
	Patients Weight	(0010,1030)	Optional
General Study	Referring Physician's Name	(0008,0090)	Optional
	Accession Number	(0008,0050)	Optional
General Series	Performing Physicians' Name	(0008,1050)	Optional
Encapsulated Document	Acquisition DateTime	(0008,002A)	
SOP Common	SOP Class UID	(0008,0016)	"1.2.840.10008.5.1.4.1.1.104.1"
	Specific Character Set	(0008,0005)	Must be "ISO_IR 100"

2.2.1.3.8 SOP Specific Conformance for 12-lead ECG Waveform Storage

HeartStation accepts a 12-lead ECG Waveform Storage object and uses the following DICOM attributes.

Table 2.2-10: 12-lead ECG Waveform Storage Attributes

Module	Attribute Name	Tag	Notes
Patient	Patient's Name	(0010,0010)	
	Patient ID	(0010,0020)	
	Patient's Birth Date	(0010,0030)	Optional
	Patient's Sex	(0010,0040)	Optional
Patient Demographics	Patient's Height	(0010,1020)	Optional
	Patients Weight	(0010,1030)	Optional
	Patient's Age	(0010,1010)	Optional
General Study	Referring Physician's Name	(0008,0090)	Optional
	Accession Number	(0008,0050)	Optional
General Series	Performing Physicians' Name	(0008,1050)	Optional
	Operators' Name	(0008,1070)	Optional
General Equipment	Manufacturer	(0008,0070)	Optional
	Institution Address	(0008,0081)	Optional
	Institution Name	(0008,0080)	Optional

Module	Attribute Name	Tag	Notes
	Station Name	(0008,1010)	Optional
	Manufacturer's Model Name	(0008,1090)	Optional
	Device Serial Number	(0018,1000)	Optional
	Software Versions	(0018,1020)	Optional
SOP Common	SOP Class UID	(0008,0016)	"1.2.840.10008.5.1.4.1.1.9.1.1"
	Specific Character Set	(0008,0005)	Must be "ISO_IR 100"
Waveform Identification	Acquisition DateTime	(0008,002A)	
Waveform Module	Waveform Sequence	(5400,0100)	
	>Waveform Originality	(003A,0004)	
	>Number of Waveform Channels	(003A,0005)	
	>Number of Waveform Samples	(003A,0010)	
	>Sampling Frequency	(003A,001A)	
	>Waveform Data	(5400,1010)	
	>Channel Definition Sequence	(003A,0200)	
	>>Channel Sensitivity	(003A,0210)	
	>>Channel Sensitivity Units Sequence	(003A,0211)	
	>>>Code Value	(0008,0010)	Must be on of mV, nV, or V
	>>Filter Low Frequency	(003A,0220)	
	>>Filter High Frequency	(003A,0221)	
	>>Notch Filter Frequency	(003A,0222)	
	>>Channel Baseline	(003A,0213)	
Waveform Annotation	Waveform Annotation Sequence	(0040,B020)	
	Measurements are obtained from the following sequence structure		
	>Concept Name Code Sequence	(0040,A043)	On sequence per measurement
	>>Code Meaning	(0008,0104)	P Duration, P Axis, PP Interval, PR Interval, QRS Axis, QRS Duration, QTc Interval, QT Interval, RR Interval, T Axis
	>Numeric Value	(0040,A30A)	Value of measurement
	Statements are obtained from the following sequence structure		
	>Concept Name Code Sequence	(0040,A043)	On sequence per statement, up to 10 supported
	>Unformatted Text Value	(0070,0006)	Statement text

2.3 Network Interfaces

HeartStation provides DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3.8 of the DICOM Standard. HeartStation inherits its TCP/IP stack from the computer system upon which it executes.

2.3.1 Physical Medium Support

HeartStation is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it is being executed.

2.4 Configuration

HeartStation configuration parameters that affect its DICOM conformance are described in this section.

2.4.1 Configuration Parameters

The following table lists all of the DICOM parameters that can be configured via the HeartStation.

Table 2.4-1: Configuration Parameter Table

Parameter	Configurable (Yes/No)	Default Value
General Parameters		
Listening Port	Yes	10104
Maximum number of simultaneous Associations	Yes	10
Time-out waiting for A-ASSOCIATE RQ on open TCP/IP connection (ARTIM timeout)	Yes	5s
Time-out waiting on an open association for the next message (DIMSE timeout)	Yes	No timeout
Time-out waiting for acceptance or rejection Response to an Association Open Request. (Application Level timeout)	No	No timeout
Time-out waiting on an open association for the next message after sending A-RELEASE RSP or A-ABORT RQ (Closing timeout)	Yes	500ms
Maximum PDU size the AE can receive	Yes	16352
Maximum PDU size the AE can send	No	65535
Pack Command and Data PDVs in one PDU	Yes	False
Support for the Basic TLS Secure Transport Connection Profile	Yes	Off
Accepted TLS Ciphers	Yes	-
DICOM Storage SCU AE Parameters		
Accepted Called AETs	Yes	Any
Accepted Calling AET	Yes	HEARTSTATION
MWL SCU AE Parameters		
Accepted Called AE Title	Yes	MWLQuery
MWL SCP IP address	Yes	0.0.0.0
DICOM MWL Port	Yes	3320
Accepted Calling AE Title	Yes	SERVERNAME
Supported SOP Class	No	Modality Worklist Information Model – FIND
Transfer Syntax	Yes	ExplicitVRLittleEndian, ImplicitVRLittleEndian

3 SUPPORT FOR EXTENDED CHARACTER SETS

HeartStation supports ISO_IR 100 (ISO 8859-1 Latin 1) as an extended character set.

4 SECURITY

4.1 Security Profiles

HeartStation supports secure DICOM communication in conformance with the Basic TLS Secure Transport Connection Profile. At default configuration, the TLS option is deactivated.

4.2 Association Level Security

HeartStation can be configured to accept Association Requests from only a limited list of Calling AE Titles.

4.3 Application Level Security

HeartStation web module can be configured to require user authentication in order to access the user interface functionalities.