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- 1 EuroSafe Poster 2014: Towards safer imaging in neonatal & paediatric radiology
- 2 P. Leblans, D. Vandenbroucke, P. Willems, "Storage Phosphors for Medical", Materials, 2011, 4, 1034-1086.
- 3 P. Leblans, L. Struye, P. Willems, "A new needle-crystalline computed radiography detector", J. Digital Imaging 13, 2012, 117-120.
- 4 EP1359204 A1 "Needle-shaped cylindrical storage phosphor crystals" nit
- 5 "Next generation MUSICA: more from each image" Interview with Piet Vuylsteke, PhD, Senior Researcher and Jan Leeuws, Business U
Manager Digital Radiography e
- 6 IEC62220-1:2003, "Medical electrical equipment - Characteristics of digital x-ray imaging devices - Part 1: Determination of th
detective quantum efficiency".
- 7 D. Vandenbroucke, P. Leblans, "CR Mammography: Image Quality Measurement and Model Calculation for Needle vs. Powder Imaging
Plate", Proceedings of the 10th International Workshop, IWDM 2010, Girona, Catalonia, Spain, June 16-18, 2010.
- 8 N. Marshall, K. Lemmens, H. Bosmans, "Physical evaluation of a needle photostimulable phosphor based CR mammography system",
Med. Phys. 39(2), 2012, 811-824
- 9 R. Schaetzing, "Management of pediatric radiation dose using Agfa computed radiography", Pediatr. Radiol. 34(3), 2004, S207-S214.
- 10 M.Cohen, D. Corea, M. Wanner, B. Karmazyn, R. Gunderman, K. Applegate, S. Jennings, Academic Radiology, 18(2), 2011, 197-198.
- 11 Gammex 610 Neonatal Phantom Users Manual
- 12 International Standard IEC 62494-1 (2008) Medical electrical equipment—exposure index of digital X-ray imaging systems
- 13 Steven Don, & Bruce Whiting & Lois Rutz & Bruce Apgar, "New digital radiography standards simplified for radiologists and
technologists," AJR:199, December 2012
- 14 Matthew Cooper, MD, Mervyn Cohen, MD; Kelly Piersall, RT; Bruce Apgar, BS (2011) "Using the exposure index to monitor radiation
exposure for portable chest radiographs in neonates". Pediatr Radiol 41:592–601
- 15 Mervyn D. Cohen & Richard Markowitz & Jeanne Hill & Walter Huda & Paul Babyn & Bruce Apgar "Quality assurance: a comparison
study of radiographic exposure for neonatal chest radiographs at 4 academic Hospitals". Pediatr Radiol 2011 Nov 6
- 16 Customer case Loma Linda, Zwanger-Pesiri radiology

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© 2014年Agfa HealthCare NV版权所有