G334(i) FIXER

G334(i) is Agfa's standard two-part universal fixer. G334(i) emits less process-specific odors while maintaining its high quality.

- Designed for high-end automatic processors
- Less emission of process-specific odors
- Ecological packaging: no separate dismissal of empty bottles
- Screw caps with built-in cutting edges for easy opening
- Bottles fit for automatic mixing units

Standard universal fixer

Agfa's two-part fixer, G334(i), excels by a lower emission of process-specific odors. This has been achieved by reducing the emission of both acetic acid and sulphur dioxide. The power and the accuracy of the fixer have not diminished and you will not, therefore, need to increase the replenishment rates of the processor.

Flexible and universal

G334(i) is suitable for use in automatic processors and will allow you to process a wide variety of X-ray films such as CURIX, SCOPIX and MAMORAY films. G334(i) fixer has been designed to be used in combination with G138i or G139 developer.

Reduced emission of chemical substances

The majority of chemical substances used in common developers and non-hardening fixers dissolve readily in water and do not excessively evaporate into the air.

G334(i) fixer, however, contains an acetic acid pH buffer and sodium sulphite to prevent the fixing salt in the fixer from decomposing. The emission of both acetic acid and sulphur dioxide using G334(i) has been greatly reduced.

E.O.S. technology compatible

G334(i) is also suitable for on-line silver recovery. In combination with the Ecologically Optimized System (E.O.S.) it is possible to obtain a minimal replenishment rate







Ecologically improved packaging

The G334(i) fixer is supplied in shipping cartons manufactured using 60 % recycled cardboard. The cartons have been provided with two handles for easy tilting and transportation. Each carton contains two two-part kits, each sufficient for preparing 25 liters (G334i) or 20 liters (G334) of ready-to-use solution. Parts A and B are packed in translucent polyethylene (PE) bottles fitted with PE ECO seals.

The screw cap and the stripe at the top of the label are color-coded for identification purposes (blue for fixer). The colored components, however, have been reduced as much as possible and the bottles show no aluminium remnants after use. Both features make the bottles easily PE-recyclable.

Easy bottle handling

You can easily open the bottles using the built-in cutting edges of the cap. Simply open the cap, place it upside down on the bottle, apply a little pressure and turn the cap clockwise to break the ECO seal.

Easy automatic mixing

You can easily mix the parts using the automatic mixing units. Remove the cap, invert the bottle with the ECO seal over the seal opener and press. The seal is automatically broken and the contents can flow into the tank on the mixing unit. Whatever the mixing procedure, automatic or manual, mixing will never result in sludge formation.

Technical specifications

GENERAL

- Suited for: all kinds of processors, except tabletop processors
- Not suited for: manual processing
- Processing cycles: all processing cycles
- Replenishment rate: ranging from 300 to 800 ml/m²

Conditions

 Make sure the processing room is well ventilated: air ventilation of 10 times the room volume per hour will prevent build-up of chemicals in the air

Packaging

- One shipping carton G334 contains: concentrate for 2 x 20 l of fixer
- One shipping carton G334i contains: concentrate for 2 x 25 l of fixer

SAFETY

 Consult the Safety Data Sheet (SDS) for information relating to safety, health and environmental aspects. Please follow the safety instructions mentioned on the bottles carefully.

For more information on Agfa, please visit our website on www.agfa.com ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2018 Agfa NV All rights reserved Published by Agfa NV Septestraat 27 - 2640 Mortsel Belgium

NE9XE GB 00201803

