G138i DEVELOPER

G138i is Agfa's standard three-part universal developer with improved consistency and stability.

- Designed to be used for automatic HT, IP or RP film processing
- Optimum results when combined with G334(i) fixer
- Ecological packaging: no separate dismissal of empty bottles

G138i for automatic processing

G138i is a three-part developer that is Agfa's standard universal developer for automatic processing in conventional systems.

G138i yields both outstanding sensitometric and physical properties. Depending on the film type used and the throughput —average processed square metres per day— replenishment may vary from 300 to 600 ml/m².

For a wide range of use

G138i can be used in all current automatic processors. HT, IP or RP processing, G138i always performs excellently. In order to obtain constant, reliable results, it is advised to use G138i developer in combination with G334(i) fixer.

G138Si starter for a perfect start

It is recommended to use G138Si starter each time you replace developer. Adding 25 ml of G138Si starting solution per liter of replenisher will allow you to obtain images featuring optimum sensitometric characteristics. Moreover, the difference in image quality during start-up time will be minimized in that way.







Screw caps with built-in cutting edges for easy opening.



Bottles fit for automatic mixing units.

Ecologically improved packaging

G138i is shipped in cartons made of 60% recycled cardboard. Two handles at both sides of the carton facilitate handling and transport.

The three parts are packed in translucent polyethylene (PE) bottles with PE ECO seals. The screw caps and the stripe at the top of the label are colored for identification purposes. Developer bottles have a red cap, except for part C, which has a transparent cap to distinguish it from part B.

The colored components, however, have been reduced as much as possible and the bottles show no aluminium remnants after use. Both aspects make the bottles easily PE-recyclable. Each carton contains concentrates for $2 \times 20 \, \text{l}$ of developer.

Easy bottle handling

You can easily open the bottles using the built-in cutting edges of the cap. Simply open the cap, place it upside down on the bottle, apply a little pressure and turn the cap clockwise to break the ECO seal.

Easy automatic mixing

You can easily mix the parts using the automatic mixing units. Remove the cap, invert the bottle with the ECO seal over the seal opener and press. The seal is automatically broken and the contents can flow into the tank on the mixing unit.

Whatever the mixing procedure, automatic or manual, mixing will never result in sludge formation.

Technical specifications

GENERAL

- Temperature adjustment of the developer by 1 °C equals speed adjustment by 2 to 3 log E points
- Suited for: all kinds of processors, except tabletop processors.
- Not suited for: manual processing
- Processing cycles:
 - CURIX HT-330 U, CURIX Compact Plus, Classic E.O.S. and Compact E.O.S.: HT, IP, RP processing cycles at 32 °C to 38 °C
 - CURIX HT-530 E.O.S.: 45 s processing cycle at 36 °C to 38 °C, 90 s and 180 s processing cycles at 32 °C to 35 °C
- Replenishment rate: ranging from 300 to 600 ml/m²
- Recommended fixer:
 - G334
 - G334i

Conditions

 Make sure the processing room is well ventilated: air ventilation of 10 times the room volume per hour will prevent build-up of chemicals in the air

Packaging

 One shipping carton G138i contains: concentrate for 2x20 liter of developer

SAFETY

Consult the Safety Data Sheet (SDS) for information relating to safety, health and environmental aspects. Please follow the safety instructions mentioned on the bottles carefully.

For more information on Agfa, please visit our website on www.agfa.com

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. DX-G and DX-M are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2018 Agfa NV All rights reserved Published by Agfa NV Septestraat 27 - 2640 Mortsel Belgium

NEJ8Y GB 00201807

