

# **AGFA HEALTHCARE**

## **IHE Integration Statement**

→ **NX 2.0.8400 / 3.0.8400**

Status: Released  
Document No. 001340  
Revision: 1.1  
NodeID LiveLink: 31341528

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfa.com/en/he/support/support_service/index.jsp">http://www.agfa.com/en/he/support/support_service/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © October, 10**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and frontpage, adjust links.

IHE Integration Statement NX 2.0.8400/3.0.8400		
Revision Number	Date	Reason for Change
1.0	September 2010	Initial version
1.1	October 2010	Add support for actor Evidence Creator in the CPI Profile

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	NX	2.0.8400 / 3.0.8400	October, 2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management	
Patient Info Reconciliation	Acquisition Modality	none	
Consistent Presentation of Images	Acquisition Modality	none	
	Evidence Creator	none	
	Print Composer	none	
Audit Trail and Node Authentication	Secure Application	none	
Consistent Time	Time Client	none	
Mammo Integration Profile	Acquisition Modality	none	
Portable Data for Imaging	Portable Media Creator	none	
Internet address for vendor's IHE information: <a href="http://www.agfa.com/healthcare/ihe">www.agfa.com/healthcare/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	<a href="http://www.agfa.com/healthcare/hl7">www.agfa.com/healthcare/hl7</a>		
DICOM	<a href="http://www.agfa.com/healthcare/dicom">www.agfa.com/healthcare/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net/">http://www.ihe-europe.net/</a>	
In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>			



## Details as of PDF Creation Date

### Document Metadata

<b>Title:</b>	001340_NX 2.0.8400/3.0.8400 IHE Integration Statement
<b>Livelink ID:</b>	31341528
<b>Version#:</b>	3
<b>Version Date:</b>	2010/10/05 05:11 PM CET
<b>Status:</b>	Approved on 2010/10/12 09:10 AM CET
<b>Owner:</b>	Peter Buytaert (awabr)
<b>Created By:</b>	Peter Buytaert (awabr)
<b>Created Date:</b>	2010/09/27 10:52 AM CET
<b>PDF Creation Date:</b>	2010/10/12 09:10 AM CET

**This document was approved by:**

### Signatures:

1. John Gibbs (ambcj) on 2010/10/12 08:25 AM CET
2. Bruno Laffin (nawyv) on 2010/10/05 08:03 PM CET
3. Koen Fabri (amlse) on 2010/10/06 09:17 AM CET

### Detailed Approver History:

- **Approval Workflow started on 2010/10/05 05:13 PM CET**
  - Approval task originally assigned to and completed by Koen Fabri (amlse) on 2010/10/06 09:17 AM CET
  - Approval task originally assigned to and completed by John Gibbs (ambcj) on 2010/10/12 08:25 AM CET
  - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2010/10/05 08:03 PM CET

### Version & Status History

Version#	Date Created	Status
3	2010/10/05 05:11 PM CET	Approved - 2010/10/12
2	2010/09/27 10:57 AM CET	Rejected - 2010/10/04
1	2010/09/27 10:52 AM CET	