

# **AGFA HEALTHCARE**

## **IHE Integration Statement**

→ **NX 2.0.8900 / 3.0.8900**

Document No. 001493  
Revision: 1.0  
NodeID Livelink: 47548470

**When printed, this is NOT a controlled copy**

## Document Information

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfahealthcare.com/global/en/main/contact/index.jsp">http://www.agfahealthcare.com/global/en/main/contact/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © 2014**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front-page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement NX 2.0.8900/3.0.8900		
Revision Number	Date	Reason for Change
1.0	November 2014	Initial version

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	NX	2.0.8900 / 3.0.8900	December, 2014
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
SWF - Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management	
PIR - Patient Info Reconciliation	Acquisition Modality	none	
ED - Evidence Documents	Acquisition Modality	none	
CPI - Consistent Presentation of Images	Acquisition Modality	none	
	Evidence Creator	none	
	Print Composer	none	
ATNA - Audit Trail and Node Authentication	Secure Application	none	
CT - Consistent Time	Time Client	none	
MAMMO - Mammo Integration Profile	Acquisition Modality	none	
PDI - Portable Data for Imaging	Portable Media Creator	none	
REM - Radiation Exposure Monitoring	Acquisition Modality	none	
Internet address for vendor's IHE information: <a href="http://www.agfahealthcare.com/ihe">www.agfahealthcare.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	<a href="http://www.agfahealthcare.com/hl7">www.agfahealthcare.com/hl7</a>		
DICOM	<a href="http://www.agfahealthcare.com/dicom">www.agfahealthcare.com/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org/">www.ihe-j.org/</a>