


AGFA HEALTHCARE

IHE Integration Statement



Drystar 5300

Status: Released
Document No.: 000955
Revision: 2.1
NodeID Livelink: 9764839

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © April, 10
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.

IHE Integration Statement Drystar 5300		
Revision Number	Date	Reason for Change
1.0	09.04.2004	Initial revision
1.1	16.12.2005	Update according new template
2.0	21.03.2007	Update according to template NodeID: 8818216 Added table 2 with latest released software version 2.0.0 C1
2.1	31.03.2010	Update Integration Profiles for SW version 2.0.0 C1

Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5300	1.82 C1	10.05.2004
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server	None	
Basic Security	Secure Node	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5300	2.0.0 C1	09.02.2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server (PS)	None	
Audit Trail and Node Authentication (ATNA)	Secure Node (SN)	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: http://www.ihe-europe.net/	
In Japan: www.jira-net.or.jp/ihe-j			



Details as of PDF Creation Date

Document Metadata

Title:	000955_Drystar 5300 IHE Integration Statement
Livelihood ID:	9764839
Version#:	6
Version Date:	2010/03/31 12:00 PM CET
Status:	Approved on 2010/04/12 11:16 AM CET
Owner:	Peter Buytaert (awabr)
Created By:	Linda De Wolf (agldw (Delete) 7337754)
Created Date:	2005/12/16 07:29 AM CET
PDF Creation Date:	2010/04/12 11:17 AM CET

This document was approved by:

Signatures:

1. Bruno Laffin (nawyv) on 2010/04/07 02:57 PM CET

Detailed Approver History:

- **Approval Workflow started on 2010/03/31 12:02 PM CET**
 - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2010/04/07 02:57 PM CET

Version & Status History

Version#	Date Created	Status
6	2010/03/31 12:00 PM CET	Approved - 2010/04/12
5	2007/03/28 11:49 AM CET	Obsolescence Check Event Completed - 2008/03/27 Approved - 2007/03/28
4	2007/03/21 01:46 PM CET	Approved - 2007/03/28
3	2007/02/27 01:21 PM CET	Approved - 2007/02/27

2	2005/12/16 07:32 AM CET	Approved - 2007/02/27
1	2005/12/16 07:29 AM CET	