

AGFA NV

IHE Integration Statement

→ **NX 3.0.22.00 / 4.0.22.00**

Document No. 001611
Revision: 3
NodeID Livelink: 69367831

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	https://medimg.agfa.com/main/contact/

Issued by:
Agfa NV
V&V Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 8219
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2019

Agfa NV

All rights reserved

REVISION RECORD

For detailed version history and version numbers, see Livelihood.

IHE Integration Statement, Global Template, NodeID: 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelihood NodeID in header and front-page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement NX 3.0.22.00 / 4.0.22.00		
Revision Number	Date	Reason for Change
1	2019-06-24	Initial version
2	2019-06-24	Add Livelihood NodeID
3	2019-06-28	Added Document Number

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA	NX	3.0.22.00 / 4.0.22.00	June, 2019
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ITI – IT Infrastructure Domain			
ATNA - Audit Trail and Node Authentication	Secure Application	none	
CT - Consistent Time	Time Client	none	
RAD - Radiology Domain			
CPI - Consistent Presentation of Images	Acquisition Modality	none	
	Evidence Creator	none	
	Print Composer	none	
ED - Evidence Documents	Acquisition Modality	none	
MAMMO - Mammo Integration Profile	Acquisition Modality	none	
PDI - Portable Data for Imaging	Portable Media Creator	none	
PIR - Patient Info Reconciliation	Acquisition Modality	none	
REM - Radiation Exposure Monitoring	Acquisition Modality	none	
SWF - Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management	
Internet address for vendor's IHE information:			
https://medimg.agfa.com/main/miscellaneous/interoperability/ihe_integration			
Links to Standards Conformance Statements for the Implementation			
HL7	https://medimg.agfa.com/main/miscellaneous/interoperability/hl7_conf_statements		
DICOM	https://medimg.agfa.com/main/miscellaneous/interoperability/dicom_connectivity		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/



Details as of PDF Creation Date

Document Metadata

Title:	001611_NX 3.0.22.00/4.0.22.00 IHE Integration Statement
Livelink ID:	69367831
Version#:	3
Version Date:	2019-06-28 11:18 AM CET
Status:	Approved on 2019-06-28 11:57 AM CET
Owner:	Dragan Jichici (Admin) (awtse_Admin)
Created By:	Koen Fabri (amlse)
Created Date:	2019-06-24 06:16 PM CET
PDF Creation Date:	2019-06-28 11:57 AM CET

This document was approved by:

Signatures:

1. Koen Fabri (amlse) on 2019-06-28 11:44 AM CET

Detailed Approver History:

- **Approval Workflow started on 2019-06-28 11:40 AM CET**
 - Approval task originally assigned to and completed by Koen Fabri (amlse) on 2019-06-28 11:44 AM CET

Version & Status History

Version#	Date Created	Status
3	2019-06-28 11:18 AM CET	Approved - 2019-06-28
2	2019-06-24 06:19 PM CET	Approved - 2019-06-25 Reviewed - 2019-06-25
1	2019-06-24 06:16 PM CET	