HE/001644

Document number 001644 - Revision 7 NodelD Livelink: 78094751

Agfa HealthCare

Page 1 of 6

# AGFA HEALTHCARE IHE Integration Statement

# → Agfa HealthCare Enterprise Imaging 8.2.x

Document No.: 001644

Revision: 7

NodeID Livelink: 78094751

When printed, this is NOT a controlled copy



Livelink ID: 78094751 Version: 7 Document Status: Approved

HE/001644 Page 2 of 6 Document number 001644 - Revision 7

NodelD Livelink: 78094751

Agfa HealthCare

# **Document Information**

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://global.agfahealthcare.com/main/contact/

Issued by: Agfa HealthCare **V&V** Connectivity Septestraat 27 B-2640 Mortsel

email: connectivity@agfa.com

Agfa HealthCare shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2022 Agfa HealthCare All rights reserved



Livelink ID: 78094751 Version: 7

HE/001644 Page 3 of 6

Document number 001644 - Revision 7

NodelD Livelink: 78094751

Agfa HealthCare

## **REVISION RECORD**

IHE Integration Statement Template NodeID : 8818216			
Revision Number	Date Reason for Change		
1.0	28 Sep 2005	Initial revision	
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.	
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.	
1.3	25 Jan 2011	Update layout + links	

IHE Integration Statement Enterprise Imaging 8.2.x			
Revision Number	Date	Reason for Change	
2-3	15-22 October 2021	Initial revision and corrections made in review cycle	
4-5	1-6 July 2022	Add TCE footnote	
6-7	1-7 december 2022	Added support for the following profiles / actors:	
		SWF.b as Image Display	
		DIFF, PERF, CXCAD, SMI as Image Manager / Archive	
		XDS-I.b as combined Imaging Document Source / Repository	
		IOCM as Image Display	
		ATNA Radiology audit trail option for Secure Application	



NodelD Livelink: 78094751

Page 4 of 6

HE/001644 Document number 001644 - Revision 7

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Enterprise Imaging <sup>1</sup>	8.2.0.x	December 2021
		8.2.1.x	July 2022
		8.2.2.x	February 2023

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:				
Integration Profiles Implemented	Actors Implemented	Options Implemented		
CARD -	CARD – Cardiology Domain			
CATH – Cardiac Catheterization Workflow	Image Display	none		
ECHO - Echocardiography Workflow	Image Display	none		
ITI – IT In	Ifrastructure Domain			
ATNA - Audit Trail and Node Authentication	Audit Record Repository	none		
ATNA - Addit Trail and Node Additionication	Secure Application	Radiology Audit Trail		
CT - Consistent Time	Time Client	none		
<b>PIX</b> / PIXV3 - Patient Identifier Cross-referencing <sup>2</sup>	Patient Identifier Cross-reference Consumer	none		
PDQ / PDQV3 – Patient Demographics Query <sup>2</sup>	Patient Demographics Consumer	none		
<b>XDS.b</b> - Cross-Enterprise Document Sharing <sup>2</sup>	Document Consumer	none		
<b>XDS-SD</b> - Cross-Enterprise Sharing of Scanned Documents <sup>2</sup>	Document Consumer	none		
RAD –	Radiology Domain			
	Image Manager / Archive	none		
ARI - Access to Radiology Information	Image Display <sup>3</sup>	Multiple Sources		
	Report Reader	none		
	Image Manager / Archive	none		
CPI - Consistent Presentation of Image	Image Display	none		
or i consistent resentation of image	Evidence Creator	none		
	Print Composer	none		
CXCAD – Chest X-Ray CAD Display	Image Manager / Archive	none		
<b>DBT</b> - Digital Breast Tomosynthesis	Image Manager / Archive	none		
-	Image Display	none		
DIFF – MR Diffusion Imaging	Image Manager / Archive	none		
ED - Evidence Document	Image Manager / Archive	none		
	Image Display	none		

<sup>&</sup>lt;sup>1</sup> Enterprise Imaging 8.2.x comprises the following components:

<sup>&</sup>lt;sup>3</sup> Some query keys not supported



Enterprise Imaging Core Server (comprised in Enterprise Imaging / Radiology Suite / Cardiology suite)
Enterprise Imaging Web Server (further referred to as Web Server for the server part and XERO Viewer for the display part)

<sup>&</sup>lt;sup>2</sup> Supported by Web Server component only

Livelink ID: 78094751 Version: 7

HE/001644 Page 5 of 6

Document number 001644 - Revision 7

NodelD Livelink: 78094751

#### Agfa HealthCare

**Evidence Creator** none FUS - Image Fusion Integration Profile Image Manager / Archive none Change Requester none **IOCM** - Imaging Object Change Management Image Manager/Archive none Image Display4 none Image Manager / Archive none Image Display KIN - Key Image Notes none **Evidence Creator** none Image Manager/Archive none MAMMO - Mammography Image Image Display<sup>5</sup> none Image Manager/Archive none NM - Nuclear Medicine Image Image Display none **Evidence Creator** none Image Display none Display none PDI - Portable Data for Imaging Portable Media Creator none Portable Media Importer none PERF - CT/MR Perfusion Imaging with Contrast Image Manager / Archive none PGP - Presentation of Grouped Procedures Image Manager / Archive none Image Manager / Archive none PIR - Patient Information Reconciliation **PPS Manager** none **REM** - Radiation Exposure Monitoring Image Manager/Archive none RWF - Reporting Workflow Image Manager/Archive none Report Reader<sup>5</sup> none **SINR** - Simple Image and Numeric Report Report Repository<sup>5</sup> none Image Manager / Archive SMI - Stereotactic Mammography Image none Image Manager / Archive none **PPS Manager** none SWF - Scheduled Workflow6 Image Display<sup>7</sup> none **Evidence Creator** none SWF.b - Scheduled Workflow.b Image Display none Export Selector8 none **TCE** – Teaching File and Clinical Trial Export Export Manager<sup>8</sup> none XCA-I - Cross-Community Access for Imaging Imaging Document Consumer<sup>2</sup> none Imaging Document Consumer<sup>2</sup> none XDS-I.b - Cross-Enterprise Document Sharing for Imaging Document Source9 none **Imaging** Document Repository9 none Internet address for vendor's IHE information: www.agfahealthcare.com/ihe

<sup>&</sup>lt;sup>9</sup> Supported as of version 8.2.2. Both actors Imaging Document Source and Document Repository are combined and cannot be dissociated



<sup>&</sup>lt;sup>4</sup> The 'configurable behavior' for handling images rejected for quality reasons is not supported

<sup>&</sup>lt;sup>5</sup> Instance level query keys not supported

<sup>&</sup>lt;sup>6</sup> SWF profile: 1 to 1 relationship between Accession Number and Requested procedure AND 1 to 1 relation between RP and SPS

<sup>&</sup>lt;sup>7</sup> Some query keys not supported

<sup>&</sup>lt;sup>8</sup> Both actors Export Selector and Export Manager are combined and cannot be dissociated

Document Status: Approved Livelink ID: 78094751 Version: 7

HE/001644 Page 6 of 6

Document number 001644 - Revision 7

NodeID Livelink: 78094751

Agfa HealthCare

Links to Standards Conformance Statements for the Implementation		
HL7	https://global.agfahealthcare.com/hl7	
DICOM	https://global.agfahealthcare.com/dicomconformance	

Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org





#### **Details as of PDF Creation Date**

#### **Document Metadata**

Title:	001644_Enterprise Imaging 8.2.x IHE Integration Statement
Livelink ID:	78094751
Version#:	7
Version Date:	2022-12-07 09:26 AM CET
Status:	Approved on 2022-12-08 09:03 AM CET
Owner:	Neil Taylor (apxzu (Delete) 65761654)
Created By:	Neil Taylor (apxzu (Delete) 65761654)
Created Date:	2021-10-15 10:25 AM CET
PDF Creation Date:	2022-12-08 09:05 AM CET

## This document was approved by:

#### Signatures:

- 1. Paul Seifert (amica) on 2022-12-07 02:22 PM CET
- 2. Bruno Laffin (nawyv) on 2022-12-08 09:02 AM CET

#### **Detailed Approver History:**

- Approval Workflow started on 2022-12-07 09:27 AM CET
  - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2022-12-08 09:02 AM CET
  - Approval task originally assigned to and completed by Paul Seifert (amica) on 2022-12-07 02:22 PM CET

### **Version & Status History**

Version#	Date Created	Status
7	2022-12-07 09:26 AM CET	Approved - 2022-12-08
6	2022-12-01 03:05 PM CET	Reviewed - 2022-12-07
5	2022-07-06 09:24 AM CET	Published - 2022-07-08 Approved - 2022-07-08
4	2022-06-30 03:27 PM CET	Reviewed - 2022-07-06

Document Status: Approved Livelink ID: 78094751 Version: 7

3		Unpublished - 2022-07-08 Published - 2021-10-26 Approved - 2021-10-26
2	2021-10-15 12:23 PM CET	Reviewed - 2021-10-22
1	2021-10-15 10:25 AM CET	

## **Applied Categories and Attributes:**

ITCo Library	
Document Type:	Conformance Statement
Category > SubCategory > Item:	Connectivity > Connectivity & Application > To Agfa
Content Manager:	Neil Taylor (apxzu (Delete) 65761654)
Summary:	Document: 001644, Revision: 5, IHE Integration Statement Enterprise Imaging 8.2.x
Language:	English
Availability:	Internet
Language Master Document:	
Software Download URL:	