



Enterprise Imaging Point-of-Care Ultrasound (POCUS)

Flexible workflows to support ultrasound
images captured at the point of care

That's life in **flow.**

AGFA 
HealthCare

Point of Care Ultrasound (POCUS) is an essential tool for ensuring your patients speedy care. AGFA HealthCare's Enterprise Imaging POCUS offers flexible interoperability and device connectivity to all DICOM ultrasound systems.

POCUS exams are performed and reported by the treating physician to promptly determine the best treatment for the patient. They often take place in fast-paced areas such as emergency departments and intensive care unit units, anesthesia or within an ambulance before the patient's arrival at the hospital. Because they are generally unscheduled, there may not

be any of the electronic order information created beforehand for traditional ultrasound exam workflows. With Enterprise Imaging POCUS, encounter workflows, multi-procedural exam acquisition, clinical and education reporting workflows are all supported by our analytical data reports, to meet your Quality Assurance needs.




“On-the-fly” exams

Storage, accessibility and quick availability are key to the effective use of ultrasound for emergency medicine. For encounter workflows, Enterprise Imaging automatically receives ADT or SUI data upon registration within the department. This allows the patient metadata to be populated automatically on the device modality worklist for quick accurate clinical workflows. Upon completion of the exam, images are stored in our centralized imaging platform and accessible via a link to our universal viewer from within your EHR.


Education and skill evaluation

As the use of ultrasound has become mainstream in emergency medicine, a need has emerged to ensure that a standard quality of care is provided. Enterprise Imaging POCUS offers an educational quality analysis module to help support your educational needs, based on the American College of Emergency Physicians (ACEP)-recommended Clinical Ultrasound Accreditation Program (CUAP). Enterprise Imaging offers residency workflows for capture, storage, reporting and overview grading of educational and clinically diagnostic exams. Analytical queries of this data are available via Enterprise Imaging's Business Intelligence module.


“With Enterprise Imaging POCUS, encounter workflows, multi-procedural exam acquisition, clinical and education reporting workflows are all supported by our analytical data reports, to meet your Quality Assurance needs.”



AGFA HealthCare Enterprise Imaging POCUS delivers:



Support for unscheduled exams, in a procedural or clinically indicated setting




Encounter workflows for creating and reporting POCUS exams ‘on the fly’



Centralized storage of all imaging data, providing easy access for the care team



Educational exam acquisition and reporting



Quality Assurance module for residency programs and competency guidelines, supported by our Business Intelligence module

Contact your AGFA HealthCare Client Executive to get started

or email enterpriseimaging@agfa.com



AGFA, the AGFA rhombus, and the AGFA HealthCare logo are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates. The Enterprise Imaging Platform logo, GRIP, GRIP Services, Imaging Health Network, Imaging Health Record, IMPAX, RUBEE, That's life in flow, XERO, and XERO Viewer are trademarks of AGFA HealthCare N.V., Belgium or its affiliates. All rights reserved. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringement. The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications that must be met by AGFA HealthCare. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. AGFA HealthCare diligently strives to provide as accurate information as possible but shall not be responsible for any typographical error.