Document number 001668 - Revision 5 NodeID Livelink: 81708848

Agfa HealthCare

Page 1 of 6

AGFA HEALTHCARE IHE Integration Statement

Agfa HealthCare Enterprise Imaging 8.3.x

Document No.: 001668

Revision: 5

NodeID Livelink: 81708848

When printed, this is NOT a controlled copy



Page 2 of 6

HE/001668

Document number 001668 - Revision 5

Agfa HealthCare

NodelD Livelink: 81708848

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://global.agfahealthcare.com/main/contact/

Issued by: Agfa HealthCare V&V Connectivity Septestraat 27 B-2640 Mortsel

email: connectivity@agfa.com

Agfa HealthCare shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2023 Agfa HealthCare All rights reserved



Livelink ID: 81708848 Version: 5

Page 3 of 6

Document number 001668 - Revision 5 NodeID Livelink: 81708848

Agfa HealthCare

HE/001668

REVISION RECORD

IHE Integration Statement Template NodeID: 8818216			
Revision Number Date Reason for Change		Reason for Change	
1.0	28 Sep 2005	Initial revision	
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.	
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.	
1.3	25 Jan 2011	Update layout + links	

IHE Integra	IHE Integration Statement Enterprise Imaging 8.3.x		
Revision Number	Date	Date Reason for Change	
2	June 2023	Initial revision for El 8.3.0.000	
3-4	July 2023	Update for El 8.3.1.000.	
		Added support of BRTO-II and MMRO-III profiles as Archive Actor	
5	December 2023	Update for EI 8.3.2.000	
		XDS-I.b Imaging Document source compatible with external Document Repository	



Document number 001668 - Revision 5 NodeID Livelink: 81708848

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Enterprise Imaging ¹	8.3.0.000	July 2023
		8.3.1.000	October 2023
		8.3.2.000	December 2023

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration

Profiles, Actors and Options listed below:				
Integration Profiles Implemented	Actors Implemented	Options Implemented		
CARD -	CARD – Cardiology Domain			
CATH – Cardiac Catheterization Workflow	Image Display	none		
ECHO – Echocardiography Workflow	Image Display	none		
ITI – IT Int	frastructure Domain			
ATNA – Audit Trail and Node Authentication	Audit Record Repository	none		
ATNA - Addit Trail and Node Addientication	Secure Application	Radiology Audit Trail		
CT – Consistent Time	Time Client	none		
PIX / PIXV3 - Patient Identifier Cross-referencing ²	Patient Identifier Cross-reference Consumer	none		
PDQ / PDQV3 – Patient Demographics Query ²	Patient Demographics Consumer	none		
XDS.b – Cross-Enterprise Document Sharing ²	Document Consumer	none		
XDS-SD – Cross-Enterprise Sharing of Scanned Documents ²	Document Consumer	none		
RAD – F	Radiology Domain			
	Image Manager / Archive	none		
ARI – Access to Radiology Information	Image Display ³	Multiple Sources		
	Report Reader	none		
	Image Manager / Archive	none		
CPI – Consistent Presentation of Image	Image Display	none		
On interest in the sentation of image	Evidence Creator	none		
	Print Composer	none		
CXCAD – Chest X-Ray CAD Display	Image Manager / Archive	none		
DBT – Digital Breast Tomosynthesis	Image Manager / Archive	none		
	Image Display	none		
DIFF – MR Diffusion Imaging	Image Manager / Archive	none		
ED – Evidence Document	Image Manager / Archive	none		
	Image Display	none		

¹ Enterprise Imaging 8.3.x comprises the following components:

³ Some query keys not supported



Page 4 of 6

Enterprise Imaging Core Server (comprised in Enterprise Imaging / Radiology Suite / Cardiology suite / VNA)
Enterprise Imaging Web Server (further referred to as Web Server for the server part and XERO Viewer for the display part)

² Supported by Web Server component only

Page 5 of 6

Document number 001668 - Revision 5 NodelD Livelink: 81708848

Agfa HealthCare

HE/001668

	Evidence Creator	nono
FUS – Image Fusion Integration Profile	Image Manager / Archive	none
- Illiage Fusion integration Frome	Change Requester	none
IOCM – Imaging Object Change Management	Image Manager/Archive	
DCM – imaging Object Change Management	Image Display ⁴	none
	Image Display Archive	none
KIN – Key Image Notes	Image Display	none
Rin – Rey illiage Notes	Evidence Creator	none
	Image Manager/Archive	none
MAMMO – Mammography Image		none
	Image Display ⁵	none
	Image Manager/Archive	none
NM – Nuclear Medicine Image	Image Display	none
	Evidence Creator	none
	Image Display	none
PDI – Portable Data for Imaging	Display	none
- Fortable Data for imaging	Portable Media Creator	none
	Portable Media Importer	none
PERF – CT/MR Perfusion Imaging with Contrast	Image Manager / Archive	none
PGP – Presentation of Grouped Procedures	Image Manager / Archive	none
DID Detient Information Decensiliation	Image Manager / Archive	none
PIR – Patient Information Reconciliation	PPS Manager	none
REM – Radiation Exposure Monitoring	Image Manager/Archive	none
RWF – Reporting Workflow	Image Manager/Archive	none
CIND Simple Image and Numeric Benert	Report Reader⁵	none
SINR – Simple Image and Numeric Report	Report Repository ⁵	none
SMI – Stereotactic Mammography Image	Image Manager / Archive	none
	Image Manager / Archive	none
SWF – Scheduled Workflow ⁶	PPS Manager	none
SWF - Scrieduled Workhows	Image Display ⁷	none
	Evidence Creator	none
SWF.b - Scheduled Workflow.b	Image Display	none
TCE Togething File and Clinical Trial Expert	Export Selector ⁸	none
TCE – Teaching File and Clinical Trial Export	Export Manager8	none
XCA-I – Cross-Community Access for Imaging	Imaging Document Consumer ²	none
What a first a second	Imaging Document Consumer ²	none
XDS-I.b – Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	none
i inaging	Document Repository ⁹	none
RO – Radiation Oncology Domain		

⁴ The 'configurable behavior' for handling images rejected for quality reasons is not supported

⁹ The Document Repository must be combined with the Imaging Document Source actor and cannot be dissociated.



⁵ Instance level query keys not supported

⁶ SWF profile: 1 to 1 relationship between Accession Number and Requested procedure AND 1 to 1 relation between RP and SPS

⁷ Some query keys not supported

⁸ Both actors Export Selector and Export Manager are combined and cannot be dissociated

HE/001668

Document number 001668 - Revision 5

NodelD Livelink: 81708848

Agfa HealthCare

Page 6 of 6

BRTO-II - Basic Radiation Therapy Objects II		Archive	none
MMRO-III - Multimodality Image Registration for Radiation Oncology 2018		Archive	none
Internet address for vendor's IHE information: https://www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7 https://www.agfahealthcare.com/hl7			
DICOM	https://www.agfahealthcare.com/dicom-conformance		

Links to general information on IHE		
IHE USA: https://www.iheusa.org	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org





Details as of PDF Creation Date

Document Metadata

Title:	001668_Enterprise Imaging 8.3.x IHE Integration Statement	
Livelink ID:	81708848	
Version#:	5	
Version Date:	2023-12-14 10:39 AM CET	
Status:	Approved on 2023-12-14 03:37 PM CET	
Owner:	Bruno Laffin (nawyv)	
Created By:	Bruno Laffin (nawyv)	
Created Date:	2023-05-30 10:30 AM CET	
PDF Creation Date:	2023-12-14 03:37 PM CET	

This document was approved by:

Signatures:

- 1. Paul Seifert (amica) on 2023-12-14 03:08 PM CET
- 2. Bruno Laffin (nawyv) on 2023-12-14 03:01 PM CET

Detailed Approver History:

- Approval Workflow started on 2023-12-14 02:59 PM CET
 - Approval task originally assigned to and completed by Paul Seifert (amica) on 2023-12-14 03:08 PM CET
 - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2023-12-14 03:01 PM CET

Version & Status History

Version#	Date Created	Status
5	2023-12-14 10:39 AM CET	Approved - 2023-12-14 Reviewed - 2023-12-14
4	2023-10-04 08:36 AM CET	Published - 2023-10-05 Approved - 2023-10-05
3	2023-09-15 04:42 PM CET	Reviewed - 2023-10-04
·		

2		Unpublished - 2023-10-05 Published - 2023-07-18 Approved - 2023-07-18 Reviewed - 2023-07-17
1	2023-05-30 10:30 AM CET	

Applied Categories and Attributes:

ITCo Library	
Document Type:	Conformance Statement
Category > SubCategory > Item:	Connectivity > Connectivity & Application > To Agfa
Content Manager:	Bruno Laffin (nawyv)
Summary:	The document 001668 is the IHE Integration Statement for Enterprise Imaging 8.3.1.000. This is the 2nd revision for the Enterprise Imaging version 8.3.x. This document will be updated for all subsequent releases of El 8.3.x until El 8.3 is release for Sale.
Language:	English
Availability:	Internet
Language Master Document:	
Software Download URL:	